

Quality Cleaning Products Since 1956

Quato 44

Disinfectant Cleaner Concentrate

Benefits:

Anionic tolerant – Removes film left by other cleaners
Non ionic detergent base – Hard water tolerant
Phosphate free – Environmentally responsible
Non staining – Will not discolour surfaces
Versatile – Cleans, disinfects and deodorizes
Safe – Non irritating at use dilutions
Registered – TPD & PMRA, EPA
Economical – Saves money
Concentrated – Dilutable

Description:

A concentrated detergent, disinfectant, sanitizer, fungicide, virucide, and odour counteractant specially formulated to meet hospital and institution needs. Effective against a broad range of both gram positive and gram-negative micro organisms and meets biocidal activity and organic soil tolerance requirements.

Caution:

- Use as directed
- Keep out of reach of children
- Keep from freezing
- Read MSDS completely before using this product
- A pre-cleaning step is required and recommended for heavily soiled surfaces

Directions & Dilutions:

Disinfection:

For hard surfaces (such as walls, floors, etc.)

Dilute 1:62 parts water (2oz/gal or 16ml/L) allow solution to wet all surfaces thoroughly and remain wet for 10 minutes. Let air-dry.

For toilet bowl:

First flush toilet, add 70ml (2.5oz) Quato 44 directly to the bowl water, swab, making sure to get under the rim, let stand 10 minutes, then flush.

Mildewstat:

To control mold and mildew on pre-cleaned hard, non-porous surfaces, dilute 1:62 parts water (2oz/gal) or 16ml/L. Wet surface completely. Let air -dry. Repeat application weekly or as required.

Sanitizing:

To sanitize pre-cleaned, rinsed, non-porous food contact surfaces, dilute 1:225 parts water, allow solution to contact surface for at least 60 seconds. Let air-dry..

Product Code & Size:

#8944-210	210 L
#8944-20	20 L
#8944-4	4 L

SWISH CANADA

1-800-461-7695
Fax: (705) 745-0220
E-mail: info@swish.ca

visit us online at

swishclean.com

24-7-365

SWISH USA

1-800-649-4149
Fax: 1-888-742-7315
E-mail: info@swish.ca